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Equipment (12) Physical Therapy Instruments and Equipment  
 General Medical Instruments: Impactor JMDN : 3285600

## MORACT Blade

### 【Warning】

1. Read these package inserts thoroughly before use.
2. These instruments are all non-sterilized. Clean and sterilize them before use.
3. When operating the instruments, use standard precautions for handling blood and other bodily fluids and follow preventive measures against infection.

### 【Contraindications and Prohibitions】

1. Do not use the instruments for any purposes other than what they are intended for.
2. Do not make any changes or modifications to the instruments.
3. Do not expose the instruments to any chemical solution.

### 【Shape, Structure and Principle】

1. Material  
Stainless steel
2. Shape, Structure

### 【Usage Purpose and Indications】

These instruments are reusable and designed to transfer impulse to target areas of the body.

### 【Operation Method or Usage Method】

Never fail to wash and sterilize the instruments after use. (It is recommended that the devices be sprayed with either diluted isopropyl alcohol, chlorhexidine or hibitane, and promptly wiped off.

### 【Precautions for Use】

1. Contact with a hard object, resulting in denting and/or deformation of the instrument tip, may reduce functionality and injure the human body during use.
2. The operator needs to have proper understanding about the patient's condition and knowledge of the correct usage techniques when using these instruments. No operators are advised to use these instruments without the appropriate skills.
3. Incorrect usage may cause bodily injuries.
4. When using these instruments, proper anatomical attention is required to effectively deliver the intended and appropriate impact on body tissue.
5. Close attention must be paid when handling the instruments to avoid unexpected and/or adverse effects.
6. Do not use any coarse abrasive agents or steel wool when cleaning the instruments.

### 【Storage Methods and Preservation】

1. Avoid high temperatures, humidity, grit and dust when storing these instruments.
2. Keep these instruments away from moisture or direct sunlight.
3. Keep these instruments out of reach of children.

### 【Maintenance and Inspection】

1. The user is responsible for usage, maintenance and inspection of these instruments.
2. Conduct inspections before and after each use, along with periodic checkups, in order to maximize instrument lifespan and safety.
3. Follow the below instructions to minimize aging of these instruments:
  - 1) Clean, sterilize and wipe off these instruments with a dry cloth immediately after use.
  - 2) If these instruments are insufficiently cleaned, do not sterilize them.
  - 3) Do not clean the instruments with any agent with high levels of acidity or alkalinity.
  - 4) Dry the instruments immediately after cleaning.

### 【Packaging】

3 units

### 【Product Complaints】

1. Contact the manufacturer (see below) for inquiries about product quality, safety, reliability, durability, efficacy and effectiveness.
2. If any of these products are responsible or thought to be responsible for having caused the death and/or serious injury of those affected and consequently receiving medical treatment, immediately contact the manufacturer by phone or fax (provided below).

### 【Marketing Authorization Holder and Manufacturer】

Marketing Authorization Holder and Manufacturer:  
 JMR co., Ltd.  
 641, Ohara, Nishikan-ku, Niigata city, Niigata  
 959-0511 Japan  
 Tel : 0256-78-7176 FAX : 0256-86-1461  
 E-mail: info@jmr-lab.com

Manufacturer: YAMAZAKI KINZOKU KOGYO Co., Ltd.

For Package Insert  
 Scan Here